

DISCOVER POWERFUL WASHING IN A SMALL FOOTPRINT



AMSCO®
2532
Single Chamber
Washer/Disinfector

Designed to support smaller CSSD's and Ambulatory Surgery Centers, the AMSCO® Single Chamber 2532 Washer/Disinfector delivers increased efficiencies – regardless of space constraints.



PRODUCTIVITY

- Increased throughput helps you meet the demands of your daily schedule
- Drive consistency with the 15 minute AutoSink cycle by providing a reliable and repeatable option for automating cleaning steps for quick turn items



PERFORMANCE

- The Full Coverage Spray Pattern uses our cross pattern spray arm nozzles to ensure complete coverage



EASE OF USE

- Capitalizing on every inch of space - the drop down loading door allows you to load and unload trays without the need of a transfer cart
- Access cycle parameters at any time using a USB port for downloading
- All racks conveniently designed for loading and unloading from all sides
- Ability to connect up to 6 lumened instruments (on a 3-level rack) per load
- Increase staff confidence with our intuitive LCD touch-screen control panel which provides clear visibility and access to critical information





GREEN SAVINGS

AMSCO 2532 Single Chamber Washer/Disinfectors are built to run efficiently with fast, effective cycles using concentrated and biodegradable chemistries

FAST, EFFECTIVE
CYCLE TIMES

MINIMIZE IMPACT ON
ENVIRONMENT

RUN
EFFICIENTLY

CONCENTRATED AND
BIODEGRADABLE
CHEMISTRIES

Size (W x H x D):

Exterior Dimensions:	28-1/2 x 78 x 28-1/2" (724 x 1981 x 724 mm)
Chamber Dimensions:	25 x 25-1/2 x 25-1/4" (635 x 648 x 642 mm)
Load Height:	29-3/4" (756 mm) above finished floor

ISO AAMI
15883
COMPLIANT

The AMSCO 2532 Washer/Disinfectors are a part of an unbroken chain of innovative solutions from STERIS

That is the power of...



For the washer/disinfector that fits your department's needs, contact your STERIS Representative or visit:
<http://www.steris.com/washing>



STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
440-354-2600 · 800-548-4873
www.steris.com