

*the* above  
curve

**V-PRO<sup>®</sup> 60**

**Low Temperature Sterilization System**

Continuing to transform the standard of low temperature sterilization, the V-PRO 60 Sterilizer delivers the flexibility and efficiency you won't find anywhere else...helping you meet the demands of today's high-performing OR's.



# Productivity

**Multiple Cycles** - Having the right cycle for the right instrument means your medical devices are ready when you need them.



**Non Lumen Cycle**



12 lb load vs. 10.7 lb load



**Flexible Cycle**



Process longer scopes,  
990mm (single & dual) vs. 850mm (single)<sup>1</sup>



**Lumen Cycle**



Process single, dual and  
triple channel devices

## Similar footprint but 2X the chamber space.



**V-PRO 60  
Sterilizer**



**Competitor's  
Chamber**

With twice the usable  
volume of the competitor's  
sterilizer, the V-PRO 60  
Sterilizer helps to  
maximize your busy OR's  
productivity.

# Ease-of-Use

The V-PRO 60 Sterilizers' easy to learn, user-friendly design can help simplify your day.



## VAPROX® HC Sterilant

59% H<sub>2</sub>O<sub>2</sub>

Same cup can be used in  
V-PRO® maX Sterilizer

## Minimal Aborts

Moisture check removes residual moisture  
from loads to save time and money



## Cycle Prompts

Makes it easy to  
match the correct  
cycle to your  
instruments



## USB Download

Gather cycle  
data instantly



## Flexible Installation

Cart or Counter  
Mounted,  
depending on your needs






# Compatibility


Plasma-free processing is gentle on devices.

## Chamber $H_2O_2$ Concentration (%)



**Lower concentration equals less potential to damage devices.**

## Chamber $H_2O_2$ Concentration<sup>2</sup>



### V-PRO Device Wizard

Quickly identify the medical devices that have been validated by STERIS for processing in the V-PRO Sterilizers



The V-PRO 60 Sterilizer  
by the Numbers:  
When Less Equals More.

**40% Less**

Consumable Cost

**24% Less**

Service Cost

---

**30% More  
Savings\***



Only STERIS provides you with an Unbroken Line of Certainty that protects your hospital from what puts it most at risk. A single, interconnected force of:

- Products
- Support
- Workflow Understanding
- Innovation
- Education

All designed to service you best when they are Connected.  
Facing unforeseen challenges together with a single vision.

This is the power of....

# A LINE OF 1

## **STERIS SecureCare<sup>SM</sup> Service Solutions**

Choose the service that fits you and your budget. With over 800 technicians worldwide, STERIS technicians are available 24/7 when you need it most.

## **CS-iQ<sup>®</sup> Sterile Processing Workflow Management Software**

includes three modules. As your needs evolve, this module can be upgraded to promote productivity, traceability and compliance across your entire department.



To learn more, contact your STERIS representative, visit [www.vpro60.com](http://www.vpro60.com) or call 1.800.548.4873.

<sup>1</sup>TERRAD NX Technical Information AD-52972 Rev A

<sup>2</sup>TERRAD NX Data Sheet AD-52368-01-4S

<sup>3</sup>Savings based on Federal Supply Schedule Price List and current as of 8/6/14.

<sup>4</sup>Over 7 year lifetime, when compared to TERRAD NX.

STERRAD 100 NX and STERRAD NX are registered trademarks of Johnson and Johnson Corporation.



STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 • USA  
440-354-2600 • 800-548-4873  
[www.steris.com](http://www.steris.com)