

Central Sterile Services

DISCOVER POWERFUL WASHING

WITH INCREASED EFFICIENCIES



One Integrated Approach to Healthcare.

 **STERIS**

Understanding the intricacies of the CSSD process has allowed **STERIS** to deliver a **holistic approach to washing:**

- Thoughtful analysis to optimize processes
- Superior cleaning chemistries to drive performance
- High-capacity water-conscious washer/disinfectors
- Verification products to help ensure positive outcomes

STERIS delivers the broadest portfolio of washer/disinfectors, all innovatively designed to increase efficiencies and ensure compliance and cost savings.


-  Productivity
-  Performance
-  Compliance
-  Ease of Use
-  Sustainability
-  Specifications



PRODUCTIVITY


AMSCO Washers/Disinfectors provide faster cycle times with low water consumption, allowing your facility to run at its highest level of performance while increasing throughput, reducing backlog and improving workflow efficiencies.

More Trays. More Cases.



Trays per 8 hour day

More Carts. More Cases.



Carts per 8 hour day



Options built to increase your productivity



Adding flexibility to meet your workflow needs, the **optional instrument cycle package** for cartwashers offers an excellent alternative for reprocessing instrument trays and loaner sets when reprocessing is at its peak.

As your needs evolve, productivity, traceability and compliance can be increased across your entire department with **CS-iQ® Sterile Processing Workflow Management Software**. Choose from three modules customized to fit your needs.



Increase productivity when washing **da Vinci® Xi® Endowrist® medical instruments** with our Vision® Single Chamber Washer.



The POWER
of AMSCO

PERFORMANCE


Full Coverage Cleaning

Our corner spray nozzles, strategically positioned on the 7053L washer/disinfector rack, deliver a full coverage clean to all instruments in their trays, eliminating the need to re-wash.



Full Coverage Spray Pattern

The ultimate cleaning technology is delivered with a cross pattern spray, corner nozzles, impingement options, and dual speed pumps. All working in concert to increase uptime and provide a more efficient process.



**STERIS 7053L washer
rack demonstrating
full chamber spray
coverage.**

**Traditional spray
technology with
minimal coverage.**

*If Trays Don't Fit
Correctly, How
Can Spray Arms
Clean Them
Properly?*



**The Belimed 290
racks don't have
enough space to
hold their trays.**

**If the tray isn't
completely on the
rack, the instruments
may not be getting
decontaminated as
intended.**



COMPLIANCE




Confidence in Dosing

Prolystica® Ultra Concentrates and HP Concentrates

are specially formulated to optimize the performance of your washer, maintain pump pressure, reduce washer strain, and provide instrument protection, thus extending the life of your investment.

Our high precision dosing system helps ensure accurate dosing, making every drop of chemistry count.

Confidence in Your Performance



Allow your department's workflow to continue as planned with a consistent, reproducible method to routinely check the performance of your automated washer/disinfector. **Verify™ All Clean Test Washer Indicators** provide you with the assurance your cycle parameters have been met.

Confidence and Certainty at the Instrument Level

VERIFY™ RESI-TEST® Cleaning Indicator uses sterile water which allows instruments to move forward to sterilization without having to be reprocessed again. This provides you with the assurance that your cleaning process was successful.



**ISO AAMI
15883
COMPLIANT**

Documentation of validated cycles serves you well during an audit.



EASE OF USE

AMSCO Washers/Disinfectors are specifically designed with you in mind. Ease of use and safety go hand-in-hand to help instill confidence in your staff which **keeps your CSSD running efficiently.**





USB port for downloading records

- Clear visibility and access to critical information provides visualization of cycle performance/status
- See cycle times from across the room with specific **color coded cycle phases** on the large display screen




Largest LED
Touch Screen
Controls on
the Market



Ease of use enables your staff to focus on critical manual reprocessing of instruments

Increase productivity with the ability to load and unload instrument trays from all sides of the rack.

Reduce the possibility of having to reprocess instrumentation and create workflow interruptions with our flexible lumen flushing ports which ensure that instruments are completely flushed during each cycle.





GREEN SAVINGS

Full Coverage Cleaning

AMSCO washers are built to run efficiently and minimize their impact on the environment. Fast, effective cycles using concentrated and biodegradable chemistries result in less water and energy consumption, and a lower total cost of ownership.

Heat
REDUCTION
Water
Waste

We develop our products to consistently deliver energy savings.


Chemical Use
SAFE
Materials
Manufacturing
Process

Promotion of
SUSTAINABLE
Designs

SAVE on UTILITIES

Vision® Single Chamber Washer uses up to

36%
LESS
water per cycle



Single Chamber Washers

Vision® 1327 Cart Washer uses up to

75%
LESS
water per cycle



Cart Washers

OUR CART WASHERS FEATURE:

- Lowest water consumption per fill of any washer on the market



- 70% recycling of water during each cycle (excludes instrument cycle)

- 10-minute cycles to reduce utility consumption compared to previous models



- Compatibility with ultra concentrated chemistries which reduces the amount of chemistry needed per cycle compared to cart washers that use non-concentrated formulations



- Recyclable packaging and biodegradable chemistry formulations



SPECIFICATIONS

Ultrasonic Cleaning Systems

	Caviwave®	Caviwave® Pro
Standard Instrument Cycle Time	4 minutes	6 minutes
Validated Automated da Vinci® Instrument Cycle	N/A	35
CaviClean Technology	N/A	Yes
Process up to 35lb tray capacity	Yes	Yes
Automated lift and lower operation	Yes	Yes
Lumen flush ports	10	10



Caviwave® Pro
Ultrasonic
Cleaning System

Single Chamber Washer/Disinfectors

	AMSCO® 2532	AMSCO® 3052	AMSCO® 5052	Vision SC®	AMSCO® 7053L
Cycle Time	30	35	25	20	27
Rack Capacity	6	10	10	10	15
Cycles Per Day	16	13	19	24	17
Trays Per Day	96	130	190	240	255
Total Water Per Cycle	21.2	47.5	38	28	39.6

Based on an 8 hour day using 5 level rack



AMSCO® 2532
Single Chamber
Washer/Disinfector



AMSCO® 3052
Single Chamber
Washer/Disinfector



AMSCO® 5052
Single Chamber
Washer/Disinfector



Vision®
Single Chamber
Washer/Disinfector

Solutions to Fit Your Needs

Making purchasing decisions while complying with the latest standards is simplified with AMSCO Washer/Disinfectors. Our wide range of sizes and functionality means AMSCO Washer/Disinfectors can easily help you meet your space, volume, and budget demands.

Multi Chamber Washer/Disinfectator

	Vision MultiChamber*
Cycle Time	12
Rack Capacity	10
Cycles Per Day	40
Trays Per Day	400
Total Water Per Cycle	38

*When running continuous workflow.

Cart Washer/Disinfectors

	Reliance 1227	Vision 1321	Vision 1327	Vision® 1300L
Cycle Time	15	10	10	10
Rack Capacity	2	1	2	3
Cycles Per Day	32	48	48	48
Trays Per Day	64	48	96	144
Total Water Per Cycle	50	14	14	14



AMSCO® 7053L
Single Chamber
Washer/Disinfectator



Vision® Multi-Chamber
Washer/Disinfectator



Vision® 1321/1327/1330L
Cart Washers/Disinfectors



1227
Cart Washer/Disinfectator

SecureCare®
Service Solutions
for Medical Equipment

Choose the service that fits you and your budget. With over 800 technicians worldwide, we are available 24/7 when you need us the most.



The Diagnostic Solutions Program diagnoses and makes timely recommendations for proper repair in the areas of water quality, chamber scale, steam components, and instrument staining.

AMSCO Washer/Disinfectors are part of an unbroken chain of innovative solutions from STERIS.

That is the power of...



da Vinci® EndoWrists®
are a registered trademark
of Intuitive Surgical Inc.

©2016 STERIS Corporation.
All rights reserved. Printed in USA.

M4015EN.2016-02 Rev. B
Printed 02/2016, XXXX

**For the Washer/Disinfecter that fits your department's needs,
contact your STERIS Representative or visit:
<http://www.steris.com/products/washing-and-decontamination-systems>**

 **STERIS**

STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
440-354-2600 · 800-548-4873
www.steris.com