



# FROM STERILIZATION TO STORAGE

POUCHCARE STERILIZATION STORAGE  
BASKETS AND DIVIDERS

A complete line of Sterilization and Storage Baskets and Dividers designed to simplify the transfer and storage of peel pouches and small packs.



# Pouchcare Sterilization and Storage Baskets and Dividers

## Flexible

- Multiple sizes both with and without baskets
- For use in Steam, Ethylene Oxide, and Vaporized Hydrogen Peroxide Processes

## Compatible

- Durable stainless steel construction
- Fits all sterilizer carts
- Washer safe

## Simple

- Quickly move pouches from sterilizer to storage
- One basket or rack for all sterilization systems
- Basket provides stackable storage when not in use

### Basket with Dividers

Part #	Description	Basket Dimensions (cm) L x W x H	Rack Dimensions (cm) L x W x H	# of Dividers
AX00072	Pouchcare Basket with Dividers-Small	14" x 11" x 4" (35 x 28 x 10)	14" x 3" x 3.5" (35 x 8 x 9)	5
AX00074	Pouchcare Basket with Dividers-Large	22" x 14" x 4" (55 x 36 x 10)	22" x 5" x 3.5" (55 x 12 x 9)	8
AX00076	Pouchcare Basket with Dividers-XL	22" x 14" x 8" (55 x 36 x 20)	22" x 5" x 7.5" (55 x 12 x 19)	8

### Baskets

Part #	Description	Quantity	Dimensions (cm) L x W x H
AX00071	Pouchcare Basket-Small	Each	14" x 11" x 4" (35 x 28 x 10)
AX00073	Pouchcare Basket-Large	Each	22" x 14" x 4" (55 x 36 x 10)
AX00075	Pouchcare Basket-XL	Each	22" x 14" x 8" (55 x 36 x 20)

### Dividers

Part #	Description	Quantity	Dimensions (cm) L x W x H	# of Dividers
AX00077	Pouchcare Dividers Only-Small	Each	13" x 5" x 3.5" (33 x 12 x 9)	5
AX00078	Pouchcare Dividers Only-Large	Each	21" x 5" x 3.5" (53 x 12 x 9)	8

Solutions designed to service you best when they are connected. Facing unforeseen challenges together with a single vision.

That is the power of...



A perfect complement to Vis-U-All® Low and High Temperature Sterilization Pouches.



\* Low and High Temperature Sterilization Pouches sold separately.

Visit [www.steris.com](http://www.steris.com) or contact your STERIS Sales Representative to learn more.