

DISCOVER POWERFUL WASHING

WITH INCREASED PRODUCTIVITY



AMSCO®
5052
Single Chamber
Washer/Disinfector

The AMSCO Single Chamber 5052 Washer/Disinfector boasts higher throughput to increase your productivity, while helping you maintain compliance.



PRODUCTIVITY

- Increased throughput with a dual wash phase allows you to focus on what's important – meeting OR demands
- Allows for flexibility when reprocessing a wide range of surgical instruments with multiple rack configurations




PERFORMANCE

- Full Coverage Spray Pattern combines our efficient cross pattern spray arms and multiple impingement options



EASE OF USE

- Increase staff confidence and safety with our intuitive LCD large touch-screen control panel which provides clear visibility and access to critical information
- Access critical cycle parameters at any time using the USB port for downloading
- Cool touch handles allow for safe and immediate handling of loads following completion of cycle
- All racks are conveniently designed for easy load and unload from all sides, minimizing hands-on time
- Innovative 180 degree reverse installation option comes standard and helps fit most space configurations





GREEN SAVINGS

AMSCO 5052 Single Chamber Washer/Disinfectors are built to run efficiently with fast, effective cycles using concentrated and biodegradable chemistries

**FAST, EFFECTIVE
CYCLE TIMES**

**MINIMIZE IMPACT ON
ENVIRONMENT**

**RUN
EFFICIENTLY**

**CONCENTRATED AND
BIODEGRADABLE
CHEMISTRIES**

Size (W x H x D):

Exterior Dimensions:	42 x 80 x 31" (1067 x 2032 x 787 mm)
Chamber Dimensions:	26-1/2 x 26-1/4 x 26" (673 x 667 x 660 mm)
Load Height:	31" (787 mm) above finished floor

**ISO AAMI
15883
COMPLIANT**

AMSCO 5052 Washer/Disinfectors are part of an unbroken chain of innovative solutions from STERIS

That is the power of...



For the washer/disinfectors that fit your department's needs, contact your STERIS Representative or visit:
<http://www.steris.com/washing>

 **STERIS**

STERIS Corporation
5960 Heisley Road
Mentor, OH 44060-1834 · USA
440-354-2600 · 800-548-4873
www.steris.com