



# The power of PRODUCTIVITY

Transforming low temperature sterilization.



## The Difference of V-PRO® Sterilizers

In the OR, ASC and the CSSD, V-PRO® Sterilizers offer the combination of increased productivity, cycle flexibility and comprehensive support you've come to expect from the leaders in reprocessing.



# Flexibility that Meets Your Demands

## Largest Chamber Capacity = Fewer Cycles to Run

	Non Lumen Cycle	Flexible Cycle	Lumen Cycle
V-PRO <sup>maX</sup> Sterilizer <i>is the most productive on the market!</i>	Up to 50 pounds (both shelves)	2 flexible endoscopes; or 1 endoscope and non-lumened load (total of 24 lbs combined)	20 stainless steel lumens [Single, dual or triple channel devices]
V-PRO <sup>60</sup> Sterilizer	Up to 12 pounds (both shelves)	1 flexible scope up to 990 mm	12 stainless steel lumens [Single, dual or triple channel devices]

## Save time and money with fewer aborted cycles.

Avoid the frustration inherent in aborted cycles. V-PRO Sterilizer's moisture-check functionality detects residual moisture in loads, thus minimizing aborted cycles and enabling greater reliability in device turnaround.



## Consistently Fast Cycle Times = Devices Ready When You Need Them



# Cost of Ownership

## V-PRO<sup>maX</sup> vs. the competition\* Sterilizer

Category	Savings by using V-PRO
Consumables (sterilant, BI, CI)	52%
Service and support	37%

**Total savings over 7 years 47%**

## V-PRO<sup>60</sup> vs. the competition\* Sterilizer

Category	Savings by using V-PRO
Consumables (sterilant, BI, CI)	40%
Service and support	24%

**Total savings over 7 years 33%**

\* As compared to U.S. government contract pricing for STERRAD 100NX and STERRAD NX in 2014







# V-PRO = Lowest Total




## Consistency Makes Operation Simple


### V-PRO Sterilizers offer you advantages that improve your bottom line:

- **GENTLER** process reduces the risk of device damage
- **FEWER** moisture related and no plasma related aborts increases productivity
- **WIDER** cycle claims and the ability to have mixed loads allows you to maximize devices per load
- **FASTER** cycle times allow for quick device turnaround

From intuitive on-screen displays to dedicated chemistries and accessories, you'll find that both operation and ongoing support were at the forefront of product design.




Large display screens feature the same standardized controls as other STERIS products, making it easy for staff to learn and operate.



**VAPROX® HC Sterilant**, an environmentally friendly and gentle hydrogen peroxide formulation, can be used in both the V-PRO maX and V-PRO 60 sterilizers. Lower H<sub>2</sub>O<sub>2</sub> concentration equals less potential to damage devices.

Our consumables are designed to make it easy to manage sterility assurance, minimize human error and protect your instruments. **Biological and chemical indicators** provide clear evidence that sterilization conditions have been achieved.





### Advantages of V-PRO Sterilizers

- Low chamber concentration makes it gentle on devices
- Self-contained cartridge limits staff exposure
- Plasma/ozone free process



The STERIS Device Testing team has the knowledge, expertise and tools to provide you confidence when reprocessing validated devices. This team works closely with medical device manufacturers to stringently test devices for material compatibility and sterilization efficacy.

The validated list of medical devices for the V-PRO Sterilizers is growing. Examples of currently validated devices include:



Endoscope	Arthroscopes Choledochoscopes Bronchoscopes	Bronchoscopes Intubation Endoscopes	Cystoscopes Laparoscopes	'da Vinci® Endoscopes Rhinolaryngoscopes	Hysteroscopes Ureteroscopes
Endoscope Accessories	Cameras	Sheaths/Adapters/Bridges	Cannulas	Fiber Optic Light Cables	Obturators
Powered Instruments	Batteries	Chuck	Drills	Saws	Handpieces
General Surgical Instruments	Curettes	Forceps	Scissors	Retractors	
Packaging	Wraps	Rigid Containers	Trays	Pouches	
Other Devices	Bougies Transducers	Defibrillation Paddles Ultrasound Probes	Laryngoscope blades and handles Gonio/Vitrectomy Lens	Flow probes	

Be notified when new device validations are complete or contact our Device Testing Team with device compatibility questions. [www.sterisdevicematrix.com](http://www.sterisdevicematrix.com)

'da Vinci endoscopes is a registered trademark of Intuitive Surgical, Inc.



CS-iQ® sterile processing workflow management software is a modular tool that helps ensure productivity, traceability.



STERIS University provides unique, quality, continuing education and training programs for both new and experienced clinicians and technicians.



The SecureCare technician network is convenient and responsive; available 24/7 to help you avoid costly downtime and keep you running smoothly. Biomed training and flexible service level options are designed to meet your needs.



# One partner. One solution.

Our philosophy involves taking a holistic approach to Central Sterile Service. Nothing is designed or engineered in a vacuum. We develop interconnected products, solutions and services that work together to refine workflow processes, increase efficiency, and provide a safe environment for patients and staff.

Our commitment to infection prevention and control is unmatched in the industry. Align with STERIS as your ONE Partner that delivers One unbroken Line of Protection.

**To learn more about V-PRO Sterilizers, contact your STERIS Representative.**

That is the power of...



STERIS Corporation  
5960 Heisley Road  
Mentor, OH 44060-1834 · USA  
800-548-4873  
[www.steris.com](http://www.steris.com)