

INFUSION RACKS

A compact, portable solution

Our selection of infusion accessories helps nursing staff to provide better care and free up vital floor space. Infusion bags and electrical syringe pumps stay out of the way yet remain within easy reach on our columns and distribution modules.

AIRport's modularity allows you to build the system you need as well as add and remove items at leisure!

• Support types:

- Rail-mounted
- Column-mounted
- Distribution module-mounted

• Support sizes:

- Single (for one infusion rack)
- Double (for two racks)

• I.V. rack configurations:

- 4 hooks (cross)
- 2 to 8 hooks (linear array)
- Basic or offset racks

• Optional height adjustment:

- With just one hand

STERIS QUALITY



• Strength:

- Supports made of epoxy-coated aluminium that resists corrosion and disinfectants
- Poles and hooks made of high-quality stainless steel



• Ergonomics:

- Single-hand height adjustment
- Offset rack possible
- Ergonomic knob for quick locking



• Design:


- Soft colours: grey and blue
- Rounded corners
- Blend in seamlessly with surgical suites and ICUs



• Examples



Rail-infusion rack clamp
Ref : TAB 741



Infusion rack
Ref : PERF 04

	Add-on rack diameter	Number of add-on racks	Length of arms		Max. range of adjustability	Rack diameter	Number of available hooks	Height adjustment	Load capacity		Weight (empty)
			Beam 1	Beam 2					Hook	Rack	
RAIL-INFUSION RACK CLAMP (25 x 10 mm)											
CLAMP 02	Ø 16	1	-	-	-	-	-	-	6 Kg	-	1 Kg
TAB 740	Ø 16	1	-	-	-	-	-	-	6 Kg	-	1 Kg
TAB 741	Ø 16	1	-	-	-	-	-	-	6 Kg	-	1,5 Kg
SUPPORT FOR INFUSION RACK FOR DISTRIBUTION COLUMN											
SUP 21	Ø 25	2	-	-	-	-	-	-	2 x 25 Kg	-	2 Kg
SUP 31	Ø 25	2	200 mm	-	200 mm	-	-	-	2 x 25 Kg	-	4,5 Kg
SUP 41	Ø 25	2	2 x 200 mm	-	2 x 200 mm	-	-	-	2 x 25 Kg	-	7 Kg
SUPPORT FOR INFUSION RACK FOR DISTRIBUTION MODULE											
SUP 12	Ø 25	1	-	-	-	-	-	-	25 Kg	-	1 Kg
SUP 22	Ø 25	1	200 mm	-	200 mm	-	-	-	25 Kg	-	3,5 Kg
SUP 32	Ø 25	1	200 mm	200 mm	400 mm	-	-	-	25 Kg	-	6 Kg
INFUSION RACK											
PERF 01	-	-	-	-	-	Ø 16	2	-	1 Kg	-	1 Kg
PERF 02	-	-	-	-	-	Ø 16	4	-	2 Kg	-	1 Kg
PERF 03	-	-	-	-	-	Ø 16	8	-	4 Kg	-	2 Kg
PERF 04	-	-	-	-	-	Ø 25	4	650 mm	2 Kg	18 Kg	2,5 Kg
PERF 05	-	-	-	-	-	Ø 25	4	650 mm	2 Kg	18 Kg	2,5 Kg
PERF 06	-	-	-	-	-	Ø 25	4	650 mm	2 Kg	18 Kg	3,5 Kg



Quad support for infusion rack for distribution column
Ref : SUP 41



Triple support for infusion rack for distribution module
Ref : SUP 32