

KEYBOARD SUPPORTS MONITOR HOLDERS

A solution for every computer interface

Our wide range of monitor holders and keyboard supports accommodates all brands of equipment and provides you with elegant and functional integrated patient management solutions.

• Monitor mounts:

- VESA (75 or 100)
- Quick-mount for Philips monitors
- Quick-mount rails (suitable for use with GE monitors)

• Range of movement:

- Combine the necessary number of arms (5 different lengths) to obtain the ideal range of movement for your workspace (192 mm to 1200 mm).

• Joints:

The number of joints determines the degrees of freedom of movement:

- Single joint
- Double joint
- Triple joint

• Optional height adjustment:

- The gas spring-assisted arm raises and lowers the monitor by up to 350 mm

• Load capacity:

- 3 to 18 kg

• Keyboard support:

- 8 models available, with or without mouse support

• Laptop support

STERIS QUALITY



• Freedom of movement:

- The monitor can be tilted sideways by 80° and up and down by 15°
- Certain models are also depth-adjustable



• Ergonomics:

- Quick mounting
- Cable management system
- Left/right sliding mouse support




• Design :

- Soft colours
- Sleek finishes
- Blends in seamlessly with surgical suites



• Examples

TECHNICAL SPECIFICATIONS



Double joint monitor holder
Ref: MON 3111




Sliding monitor holder
Ref: MON 3133



Monitor holder for Philips monitors
Ref: MON 3144



Monitor holder dual hinged joints +
keyboard support with mouse
Ref: MON 3163



Keyboard support
and mouse support
Ref: KEY 111

	Arm length		Max. range of adjustability	Monitor rotation range	Monitor tilt range	Vertical reach	Depth adjustment	Cable management system	Monitor bracket	Keyboard support	Mouse support	Load capacity	Weight (empty)
	Beam 1	Beam 2											
MONITOR HOLDER													
MON 1111	-	-	Stationary	± 80°	± 15°	-	-	-	VESA 75 / 100	-	-	8 Kg	1 Kg
MON 2111	192 mm	-	Stationary	± 80°	± 15°	-	-	Yes	VESA 75 / 100	-	-	8 Kg	2 Kg
MON 3111	320 mm	-	320 mm	± 80°	± 15°	-	-	Yes	VESA 75 / 100	-	-	8 Kg	2,25 Kg
MON 4111	320 mm	320 mm	640 mm	± 80°	± 15°	-	-	Yes	VESA 75 / 100	-	-	8 Kg	7 Kg
MON 3311	400 mm	-	400 mm	± 80°	± 15°	± 250 mm	-	-	VESA 75 / 100	-	-	8 Kg	7,5 Kg
MON 1113	-	-	Stationary	± 80°	± 15°	-	-	-	VESA 75 / 100	-	-	14 Kg	1,25 Kg
MON 2113	192 mm	-	Stationary	± 80°	± 15°	-	-	Yes	VESA 75 / 100	-	-	14 Kg	2 Kg
MON 3113	320 mm	-	320 mm	± 80°	± 15°	-	-	Yes	VESA 75 / 100	-	-	14 Kg	2,75 Kg
MON 4113	320 mm	320 mm	640 mm	± 80°	± 15°	-	-	Yes	VESA 75 / 100	-	-	14 Kg	7,25 Kg
MON 3313	400 mm	-	400 mm	± 80°	± 15°	± 250 mm	-	-	VESA 75 / 100	-	-	14 Kg	7,75 Kg
MON 2123	192 mm	-	Stationary	± 80°	± 15°	-	Yes	Yes	VESA 75 / 100	-	-	14 Kg	2,25 Kg
MON 3123	320 mm	-	320 mm	± 80°	± 15°	-	Yes	Yes	VESA 75 / 100	-	-	14 Kg	2,75 Kg
MON 4123	320 mm	320 mm	640 mm	± 80°	± 15°	-	Yes	Yes	VESA 75 / 100	-	-	14 Kg	7,25 Kg
MON 3323	400 mm	-	400 mm	± 80°	± 15°	± 250 mm	Yes	-	VESA 75 / 100	-	-	14 Kg	7,75 Kg
MON 1115	-	-	Stationary	± 80°	± 15°	-	-	-	VESA 75 / 100	-	-	18 Kg	1,5 Kg
MON 3211	750 mm	-	750 mm	± 80°	± 15°	± 350 mm	-	Yes	VESA 75 / 100	-	-	3/8 Kg	7,75 Kg
MON 3212	750 mm	-	750 mm	± 80°	± 15°	± 350 mm	-	Yes	VESA 75 / 100	-	-	10 Kg - 13 Kg	8,5 Kg
MON 4211	750 mm	450 mm	1200 mm	± 80°	± 15°	± 350 mm	-	Yes	VESA 75 / 100	-	-	3 Kg - 8 Kg	10,25 Kg
MON 4212	750 mm	450 mm	1200 mm	± 80°	± 15°	± 350 mm	-	Yes	VESA 75 / 100	-	-	10 Kg - 13 Kg	11 Kg
SLIDING MONITOR HOLDER													
MON 2133	192 mm	-	Stationary	± 80°	± 15°	-	-	Yes	Sliding	-	-	14 Kg	2 Kg
MON 3133	320 mm	-	320 mm	± 80°	± 15°	-	-	Yes	Sliding	-	-	14 Kg	2,25 Kg
MON 4133	320 mm	320 mm	640 mm	± 80°	± 15°	-	-	Yes	Sliding	-	-	14 Kg	7 Kg
MON 3333	400 mm	-	400 mm	± 80°	± 15°	± 250 mm	-	-	Sliding	-	-	14 Kg	7,5 Kg
MONITOR HOLDER FOR PHILIPS MONITORS													
MON 2144	128 mm	-	Stationary	± 80°	± 15°	-	-	-	Philips	-	-	16 Kg	2,25 Kg
MON 3144	224 mm	-	224 mm	± 80°	± 15°	-	-	-	Philips	-	-	16 Kg	2,75 Kg
MON 3344	400 mm	-	400 mm	± 80°	± 15°	± 250 mm	-	-	Philips	-	-	16 Kg	7,75 Kg
MON 2154	128 mm	-	Stationary	± 80°	± 15°	-	-	-	Philips	-	-	16 Kg	2,5 Kg
MON 3154	224 mm	-	224 mm	± 80°	± 15°	-	-	-	Philips	-	-	16 Kg	3 Kg
MON 3354	400 mm	-	400 mm	± 80°	± 15°	± 250 mm	-	-	Philips	-	-	16 Kg	8 Kg
MONITOR AND KEYBOARD HOLDER													
MON 3163	320 mm	-	320 mm	± 80°	± 15°	-	-	Yes	VESA 75 / 100	Yes	Yes	14 Kg	7,25 Kg
MON 4163	320 mm	400 mm	720 mm	± 80°	± 15°	± 250 mm	-	Yes	VESA 75 / 100	Yes	Yes	14 Kg	11,5 Kg
KEYBOARD HOLDER													
KEY 121	-	-	Stationary	-	-	-	-	-	-	Yes	Yes	3 Kg	8,5 Kg
KEY 111	468 mm	-	Stationary	-	-	-	-	Yes	-	Yes	Yes	5 Kg	3,75 Kg
KEY 311	544 mm	-	-	-	-	-	-	Yes	-	Yes	Yes	5 Kg	4,25 Kg
KEY 411	544 mm	348 mm	892 mm	-	-	-	-	Yes	-	Yes	Yes	5 Kg	8,75 Kg
KEY 131	405 mm	-	405 mm	-	-	-	-	Yes	-	Yes	Yes	5 Kg	3,75 Kg
KEY 112	-	-	-	-	-	-	-	Yes	-	Yes	-	5 Kg	1 Kg
LAPTOP HOLDER													
KEY 113	-	-	-	-	-	-	-	-	-	Yes	Yes	10 Kg	4 Kg