

## LAMPS

### A whole new way of lighting patients' bedsides

AIRport lamps are a functional and elegant solution for examining patients, reading and writing records, carrying out prep work, and performing other tasks at patients' bedsides.

They deliver high-quality light and protect patients and staff from electrical hazards. Their sleek design blends in seamlessly with patients' bedsides.

#### • Lamp types:

- Examination
- Multifunction: the same lighthead combined with a powerful examination lamp and a soft reading lamp

#### • Lighting technology:

- Halogen
- LED

#### • Arm types:

- Flexible
- Double joint
- Triple joint

#### • Brightness levels:

- From 22,500 lux at 500 mm to 48,000 lux at 500 mm

## STERIS QUALITY



#### • Safety:

- Reinforced protective front glass
- Protected transformer and wiring



#### • Ergonomics:

- Handle for easy and safe positioning
- Switch on shade
- Easy bulb replacement



#### • Design :

- Rounded corners for easy cleaning and disinfection
- Blends in seamlessly with hospital environments and ICUs



## • Examples



Halogen examination lamp  
Ref : LAMP 01



Halogen multifunction lamp  
Ref : LAMP 05

	Power consumption	Illumination intensity	Field diameter	Arm type
<b>HALOGEN EXAMINATION LAMP</b>				
LAMP 01	50 W	48 000 Lux	500 mm	Flexible
LAMP 02	50 W	48 000 Lux	500 mm	Double joint
LAMP 03	50 W	48 000 Lux	500 mm	Triple joint
<b>HALOGEN MULTIFUNCTION LAMP</b>				
LAMP 04	50 W	48 000 Lux - 1 700 Lux	500 mm - 500 mm	Flexible
LAMP 05	50 W	48 000 Lux - 1 700 Lux	500 mm - 500 mm	Double joint
LAMP 06	50 W	48 000 Lux - 1 700 Lux	500 mm - 500 mm	Triple joint
<b>LED EXAMINATION LAMP</b>				
LAMP 07	3 x 1 W	22 500 Lux	500 mm	Flexible
LAMP 08	3 x 1 W	22 500 Lux	500 mm	Double joint
LAMP 09	3 x 1 W	22 500 Lux	500 mm	Triple joint



LED examination lamp  
Ref : LAMP 09