

SHELVES DRAWERS

Build the ideal storage unit for your needs!

With 38 models of shelves and drawers to choose from, you can create the functional and stylish storage unit you need to hold your instruments, monitors, records and more.

• Sizes:

Two shelf sizes:

- **Standard model:** 46 cm wide by 36 cm deep
- **Large model:** 63 cm wide by 46 cm deep

• Drawers:

- Our shelves can accommodate up to two drawers. Each drawer has a load capacity of 3 kg.

• Mouse support:

- This ingenious addition for standard shelves slides easily to accommodate left-handed and right-handed users.

• Handles:

- Easy to grasp

• Pneumatic brake control:

- Releases each beam for easy and precise positioning

• Rails:

- For attaching many accessories, from suction bottles to waste bins. Our rails have a load capacity of 5 kg.

STERIS QUALITY



• Strength:

- Total load capacity of up to 50 kg



• Compactness:

- Maximum storage capacity with a small footprint



• Design :


- The slight inward curve along the edges prevents objects from falling and ensures easy and thorough cleaning
- The soft, restrained colours blend in seamlessly with hospital environments




• Examples




Compact shelf with rails and handles
Ref : SHE S15



Compact shelf with mouse support
and handles
Ref : SHE S23



Compact shelf with one drawer,
rails and handles
Ref : SHE S35



Compact shelf with two drawers,
rails and handles
Ref : SHE S45

TECHNICAL SPECIFICATIONS

	Overall dimensions			Platform area		Interior area			Load capacity				Weight	Options		
	With	Depth	Height	Shelf		Drawer			Shelf + Rails	Rail	1 st Drawer	2 nd Drawer		Rail	Handle	Control Pneumatic Brake
				With	Depth	With	Depth	Height								
SHELF																
SHE S11	450 mm	430 mm	32 mm	400 mm	350 mm	-	-	-	20 Kg	-	-	-	5 Kg	-	-	-
SHE S12	450 mm	430 mm	32 mm	400 mm	350 mm	-	-	-	20 Kg	5 Kg to 250 mm	-	-	7 Kg	2	-	-
SHE S13	450 mm	430 mm	32 mm	400 mm	350 mm	-	-	-	20 Kg	-	-	-	5,5 Kg	-	2	-
SHE S14	450 mm	430 mm	32 mm	400 mm	350 mm	-	-	-	20 Kg	-	-	-	6,5 Kg	-	2	1
SHE S15	450 mm	430 mm	32 mm	400 mm	350 mm	-	-	-	20 Kg	5 Kg to 250 mm	-	-	7,5 Kg	2	2	-
SHE S16	450 mm	430 mm	32 mm	400 mm	350 mm	-	-	-	20 Kg	5 Kg to 250 mm	-	-	8,5 Kg	2	2	1
SHE L11	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	50 Kg	-	-	-	7 Kg	-	-	-
SHE L12	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	50 Kg	5 Kg to 250 mm	-	-	9 Kg	2	-	-
SHE L13	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	50 Kg	-	-	-	7,5 Kg	-	2	-
SHE L14	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	50 Kg	-	-	-	9 Kg	-	2	1
SHE L15	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	50 Kg	5 Kg to 250 mm	-	-	9,5 Kg	2	2	-
SHE L16	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	50 Kg	5 Kg to 250 mm	-	-	11 Kg	2	2	1
SHELF WITH MOUSE PAD																
SHE S21	450 mm	430 mm	32 mm	400 mm	350 mm	-	-	-	20 Kg	-	-	-	6 Kg	-	-	-
SHE S23	450 mm	430 mm	32 mm	400 mm	350 mm	-	-	-	20 Kg	-	-	-	6,5 Kg	-	2	-
SHELF WITH 1 DRAWER																
SHE S31	450 mm	430 mm	128 mm	400 mm	350 mm	340 mm	240 mm	80 mm	15 Kg	-	3 Kg	-	7 Kg	-	-	-
SHE S32	450 mm	430 mm	128 mm	400 mm	350 mm	340 mm	240 mm	80 mm	15 Kg	5 Kg to 250 mm	3 Kg	-	9 Kg	2	-	-
SHE S33	450 mm	430 mm	128 mm	400 mm	350 mm	340 mm	240 mm	80 mm	15 Kg	-	3 Kg	-	7,5 Kg	-	2	-
SHE S34	450 mm	430 mm	128 mm	400 mm	350 mm	340 mm	240 mm	80 mm	15 Kg	-	3 Kg	-	7,5 Kg	-	2	1
SHE S35	450 mm	430 mm	128 mm	400 mm	350 mm	340 mm	240 mm	80 mm	15 Kg	5 Kg to 250 mm	3 Kg	-	9,5 Kg	2	2	-
SHE S36	450 mm	430 mm	128 mm	400 mm	350 mm	340 mm	240 mm	80 mm	15 Kg	5 Kg to 250 mm	3 Kg	-	9,5 Kg	2	2	1
SHE L31	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	45 Kg	-	3 Kg	-	9 Kg	-	-	-
SHE L32	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	45 Kg	5 Kg to 250 mm	3 Kg	-	11 Kg	2	-	-
SHE L33	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	45 Kg	-	3 Kg	-	9,5 Kg	-	2	-
SHE L34	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	45 Kg	-	3 Kg	-	9,5 Kg	-	2	1
SHE L35	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	45 Kg	5 Kg to 250 mm	3 Kg	-	11,5 Kg	2	2	-
SHE L36	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	45 Kg	5 Kg to 250 mm	3 Kg	-	11,5 Kg	2	2	1
SHELF WITH 2 DRAWERS																
SHE S41	450 mm	430 mm	224 mm	400 mm	350 mm	340 mm	240 mm	80 mm	10 Kg	-	3 Kg	3 Kg	9 Kg	-	-	-
SHE S42	450 mm	430 mm	224 mm	400 mm	350 mm	340 mm	240 mm	80 mm	10 Kg	5 Kg to 250 mm	3 Kg	3 Kg	11 Kg	2	-	-
SHE S43	450 mm	430 mm	224 mm	400 mm	350 mm	340 mm	240 mm	80 mm	10 Kg	-	3 Kg	3 Kg	9,5 Kg	-	2	-
SHE S44	450 mm	430 mm	224 mm	400 mm	350 mm	340 mm	240 mm	80 mm	10 Kg	-	3 Kg	3 Kg	9,5 Kg	-	2	1
SHE S45	450 mm	430 mm	224 mm	400 mm	350 mm	340 mm	240 mm	80 mm	10 Kg	5 Kg to 250 mm	3 Kg	3 Kg	11,5 Kg	2	2	-
SHE S46	450 mm	430 mm	224 mm	400 mm	350 mm	340 mm	240 mm	80 mm	10 Kg	5 Kg to 250 mm	3 Kg	3 Kg	11,5 Kg	2	2	1
SHE L41	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	40 Kg	-	3 Kg	3 Kg	11 Kg	-	-	-
SHE L42	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	40 Kg	5 Kg to 250 mm	3 Kg	3 Kg	13 Kg	2	-	-
SHE L43	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	40 Kg	-	3 Kg	3 Kg	11,5 Kg	-	2	-
SHE L44	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	40 Kg	-	3 Kg	3 Kg	11,5 Kg	-	2	1
SHE L45	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	40 Kg	5 Kg to 250 mm	3 Kg	3 Kg	13,5 Kg	2	2	-
SHE L46	650 mm	630 mm	64 mm	600 mm	450 mm	-	-	-	40 Kg	5 Kg to 250 mm	3 Kg	3 Kg	13,5 Kg	2	2	1