

XLED™: Tailor-made Surgical Lighting

Each medical specialty requires particular and specific installations.


The modular design of XLED™ surgical lighting with multiple associated products brings a unique architectural, tailor-made solution which precisely responds to both your current and future needs.

Did You Know?

Designed and developed by STERIS, XTND® suspensions are easy to install and offer a wide range of flexible and evolutive solutions: any time after installation, add one or more additional spindles to suit your needs.

1. Choose your suspension

- Choice of height (short suspensions for low ceilings)
- Choice of arm length (long arms for hybrid rooms)
- Central or dual mounted installation, up to 5 arms



2. Choose the number of spots and light heads



1 spot



2 spots



3 spots



4 spots

From 1 to 4 spots and up to 3 light heads on a single arm
Fixed or adjustable pattern diameter



Create your Future Lighting System



By choosing the number of light heads and spots, then by varying the intensity and colour temperature, you create a light solution perfectly adapted to every situation.



3. Choose the colour temperature: fixed or adjustable



4 400 K



3 500 K - 4 000 K
4 500 K - 5 000 K

4. Add cameras and monitors



HD camera handle



Spring arm and simple or dual monitor holder



Monitors available in HD or 4K, from 31 to 65 inch